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### HSL-Abutments und HSL-Kappen

B-, C-, D-, E-, F-, H-, I-, K-, L-, N-, R-, S-, T- und Y-Serie

# CE0483

#### Indikation:

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# €0483

### Indication:

For the production of prosthetic components on implants. For casting universal, casting or bake-on alloys containing noble metals.

The castable gold abutments of the B-Series are indicated for Sky® implants, manufactured by Bredent Medical®

The current implant diameters has to be attended.

The castable gold abutments of the C-Series are indicated for Camlog® implants.

The current implant diameters has to be attended.

The castable gold abutments of the D-Series are indicated for Conelog® implants.

The current implant diameters has to be attended.

The castable gold abutments of the E-Series are indicated for Replace Select® implants, manufactured by Nobel Biocare®.

The current implant diameters has to be attended.

The castable gold abutments of the F-Series are indicated for Noble Active® implants, manufactured by Nobel Biocare®.

The current implant diameters has to be attended.

The castable gold abutments of the H-Series are indicated for Certain® implants, manufactured by Biomet 3i®.

The current implant diameters has to be attended.

The castable gold abutments of the I-Series are indicated for Osseotite® implants manufactured by Biomet 3i®

The current implant diameters has to be attended.

The castable gold abutments and castable gold caps of the K-Series are indicated for Brånemark® implants, manufactured by Nobel Biocare®

The current implant diameters has to be attended.

The castable gold abutments of the L-Series are indicated for Bone Level® implants, manufactured by Straumann®.

The current implant diameters has to be attended.

The castable gold abutments of the N-Series are indicated for SynOcta® implants, manufactured by Straumann®

The current implant diameters has to be attended.

The castable gold abutments of the R-Series are indicated for Tapered Screw-Vent® implants, manufactured by Zimmer®, MIS® implants (Internal Hexagon) and BioHorizons® implants (Internal Hexagon)®.

The current implant diameters has to be attended.

The castable gold abutments of the S-Series are indicated for OsseoSpeed® implants, manufactured by Astra Tech®

The current implant diameters has to be attended.

The castable gold abutments of the T-Series are indicated for Frialit® implants, manufactured by Dentsply-Friadent®

The current implant diameters has to be attended.

The castable gold abutments of the Y-Series are indicated for Ankylos® implants, manufactured by Dentsply-Friadent®

The current implant diameters has to be attended.

The castable gold caps of the B, C, D, E, F, H, I, K, L, N, R, S, T, Y series are compatible with the MedentiBase-Abutment of each series.

The castable gold abutments and castable gold caps are delivered with screw included.

### Composition:

Castable base:

Gold (Au) 60.00% Platinum (Pt) 19.00% Palladium (Pd) 20.00% Iridium (Ir) 1.00%

Bake-on plastic sleeve: POM

Abutment screw:

Ti6Al4V, medical grade 5, ASTM 136

Technical data: Melting interval 1400 – 1490 °C Vickers hardness HV 5/30 w 145

s 225 a 240

(w) Soft annealed 1000 °C/30 min

(s ) Self-hardening by slow cooling of the cast

(a ) Hardening/annealing 700 °C/30 min, air cooling

### Contraindication:

The castable gold abutments and castable gold caps of each series can only be combined with the corresponding compatible implant system/MedentiBASE-Abutment, e.g. the E-Series only in combination with the implant system Replace Select®

No inappropriate in connection geometry abutments should be used. Any post-processing at the connection geometry to the implant system/MedentiBASE-Abutment may result in fitting inaccuracies prohibiting further use. All castable gold abutments and castable gold caps are intended for single use. There should be no other than that of the manufacturer's recommended torque be used.

Ncm Abutment

15	Bridge screw and castable gold cap
15	F-Series 3,0; S-Series 3,0
20	C-Series D-Series H-Series
25	B-Series; S-Series 3,5/4,0 4,5/5,0; T-Series; Y-Series
30	R-Series
35	I-Series; K-Series; L-Series; N-Series; E-Series; F-Series NP/RP

### Handling:

The castable gold abutment and castable gold cap is made of a base of castable alloy and a plastic sleeve that burns out without leaving a residue. The transparent plastic sleeve can be shortened. Modelling is with standard techniques.

The minimum wall thickness of 0.4 mm must be observed.

The connection geometry and platform to the implant must be completely free from plastic, wax and grease to ensure that there is no casting in this area.

### Embedding

Only phosphate-bonded (gypsum free) investment compounds recommended for noble metal casting may be used. It is especially important to ensure that there are no air bubbles in the investment.

### Preheating

Follow the processing instructions issued by the manufacturer of the investment compound. Observe the pre-heating temperature specified by the manufacturer of the alloy and increase this by 50 °C. This preheating temperature should be held for 45 minutes, for instance, in a 3x casting cylinder to ensure complete casting.

### Melting/casting

The casting temperature of the bake-on alloy must be significantly over 1400 °C (solidus temperature of the castable gold abutments and castable gold caps) to prevent it from melting. Observe the instructions for the casting equipment.

### Cooling

To prevent the cast object from being subject to stresses, cool the casting cylinder down to room temperature.

### Divesting/Stripping

The investment compound may only be removed with max. 2 bar pressure with glass polishing beads to ensure that the dimensions of the connection geometry between the implant and the MedentiBASE abutment remain unchanged. The cooled cast object may then be placed in a warm stripping solution.

### Finishing

The cast object can be finished with ceramic bonded discs / stones or with cross cut tungsten carbide burs. We recommend screwing the castable gold abutment / castable gold cap into a laboratory implant to protect the geometry of the connection during finishing.

### Veneering

If the cast castable gold abutment / castable gold cap is to be veneered, please observe the features of the ceramic compound (CTE value) and the bonding alloy. The ceramic compound may not come into contact with the alloy of the castable gold abutment / castable gold cap, due to the different CTE values (cracks, chipping).

**Safety notice:** Metal dust is harmful to health. When machining and sandblasting use a suction with a fine dust filter that is usual in practice and wear protective goggles as well as a face mask.

### Side effects:

Allergies or sensibilities in connection with the alloy cannot be ruled out in very rare instances..

### Interactions:

Various alloy types in the same mouth cavity can lead to galvanic reactions in the event of occlusal or proximal contact.

### Notice:

The information conveyed orally, in writing or in practical seminars is based on tests and experience and can therefore only be considered to be standard values. Our products are subject to constant further development. In this connection we reserve the right to make product changes in relation to the design and composition.

### Cleaning and sterilization:

All castable gold abutments and castable gold caps B, C, D, E, F, H, I, K, L, N, R, S, T and Y series are supplied non-sterile. The castable gold abutments and castable gold caps must be prepared as follows before their first use. They must be cleaned and sterilized before every use. Visually inspect to ensure they are intact and clean. If contamination is visible on the castable gold abutments and castable gold caps, clean and sterilize again until there is no visible contamination.

### Cleaning:

Automatic cleaning procedure based on Vario TD program:

- Cold water rinse 1 min.
- Cleaning at 55°C/131°F (± 2°C/35,6°F) for 5 min. with cleaning agent neodisher® MediClean (0.5% v/v)
- Neutralization with cold water for 2 min.
- Cold water rinse 1 min.

### Sterilization:


Seal in a suitable autoclave foil. Steam sterilization with fractioned vacuum process:

132°C (270°F) for 4 minutes, drying time: 20 minutes

### Warning:

The product has to be inspected prior usage. The packing must be sealed without any visible damage. The following descriptions are not sufficient for the immediate use of the castable gold abutments and castable gold caps. Dental skills and prior instruction of how to use the castable gold abutments and castable gold caps are at any rate required.

### Signs and Symbols:

-  For single use only
-  LOT number
-  Reference number
-  Manufacturer
-  Attend instruction for use
-  Non sterile

### IFU\_M0004\_REV.F\_2013-11-06

 Medentika® GmbH  
Hammweg 8–10  
D-76549 Hügelsheim/Germany  
Tel: +49-7229-69912-0  
Fax: +49-7229-69912-20  
E-Mail: info@medentika.de  
www.medentika.de

Products indicated with ® are registered brand names of the manufacturers.