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Scanbodies

B-, C-, D-, E-, F-, H-, I-, K-, L-, N-, R-, S-, T- und Y-Serie

CE 0483

Indikation:

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CE0483

Indication:

The scanbody is used to detect the exact position of the implant in the mouth or in the model to transfer it in the correct position of a 3D model.

The scanbodies of the B-Series are indicated for Sky[®] implants, manufactured by Bredent Medical[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies of the C-Series are indicated for Camlog[®] implants.

The current implant diameters and connection geometries has to be attended.

The scanbodies of the D-Series are indicated for Conelog[®] implants.

The current implant diameters and connection geometries has to be attended.

The scanbodies of the E-Series are indicated for Replace Select[®] implants, manufactured by Nobel Biocare[®].

The current implant diameters and connection geometries has to be attended.

The scanbodies of the F-Series are indicated for Nobel Active[®] implants, manufactured by Nobel Biocare[®].

The current implant diameters and connection geometries has to be attended.

The scanbodies of the H-Series are indicated for Certain[®] implants, manufactured by Biomet 3i[®].

The current implant diameters and connection geometries has to be attended.

The scanbodies of the I-Series are indicated for Osseotite[®] implants manufactured by Biomet 3i[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies of the K-Series are indicated for Brånemark[®] implants, manufactured by Nobel Biocare[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies of the L-Series are indicated for Bone Level[®] implants, manufactured by Straumann[®].

The current implant diameters and connection geometries has to be attended.

The scanbodies of the N-Series are indicated for SynOcta[®] implants, manufactured by Straumann[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies of the R-Series are indicated for Tapered Screw-Vent[®] implants, manufactured by Zimmer[®], MIS[®] implants (Internal Hexagon) and BioHorizons[®] implants (Internal Hexagon)[®].

The current implant diameters and connection geometries has to be attended.

The titanium bases and scanbodies of the S-Series are indicated for OsseoSpeed[®] implants, manufactured by Astra Tech[®]

The current implant diameters and connection geometries has to be attended.

The titanium bases and scanbodies of the T-Series are indicated for Frialit[®] implants, manufactured by Dentsply-Friadent[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies of the Y-Series are indicated for Ankylos[®] implants, manufactured by Dentsply-Friadent[®]

The current implant diameters and connection geometries has to be attended.

The scanbodies are provided with screw included.

Composition:

Abutment screw: Ti6Al4V, medical grade 5, ASTM 136
Scanbody 1. Generation: Polyether-ether-ketone, Peek
Scanbody 2. Generation:
Ti6Al4V, medical grade 5, ASTM 136; special coated
Scanbody MedentiBASE:
Ti6Al4V, medical grade 5, ASTM 136; special coated
Scanbody Multi-unit:
Ti6Al4V, medical grade 5, ASTM 136; special coated

Contraindication:

The scanbodies of each series can only be combined with the corresponding compatible implant system, e.g. the E-Series only in combination with the implant system Replace Select[®]. No inappropriate in connection geometry scanbodies should be used. Any post-processing at the connection geometry to the implant may result in fitting inaccuracies prohibiting further use. Furthermore, the scanbodies may not be ground and their shape must not be changed, as this is essential to ensure an accurate scan result.

Side effects:

Allergies or sensitivities in connection with the alloy cannot be ruled out in very rare instances.

Interactions:

Various alloy types in the same mouth cavity can lead to galvanic reactions in the event of occlusal or proximal contact.

Notice:

The information conveyed orally, in writing or in practical seminars is based on tests and experience and can therefore only be considered to be standard values. Our products are subject to constant further development. In this connection we reserve the right to make product changes in relation to the design and composition.

Cleaning and sterilization:

The Scanbody 2. Generation, Scanbody MedentiBASE and Scanbody Multi-unit of the B, C, D, E, F, H, I, K, L, N, R, S, T, Y series are delivered in a non-sterile state in the appropriate packaging. The products has to be cleaned and sterilized according to following details:

Cleaning:

Automatic cleaning procedure based on Vario TD program:

- Cold water rinse 1 min
- Cleaning at 55°C/131°F (± 2°C/35,6°F) for 5 min with cleaning agent neodisher[®] MediClean (0.5% v/v)
- Neutralization with cold water for 2 min
- Cold water rinse 1 min

Sterilization:

Steam sterilize the scanbodies using the cycle similar to instrumentation sterilization such as 4 minutes at 132°C/270°F.

Warning:

The product has to be inspected prior usage. The packing must be sealed without any visible damage. The following descriptions are not sufficient for the immediate use of the product. Dental skills and prior instruction of how to use the product are at any rate required.

Signs and Symbols:

-  For single use only
-  LOT LOT number
-  REF Reference number
-  Manufacturer
-  Attend instruction for use
-  Non sterile

IFU_M0009_REV.A_2013-07-10

Medentika[®] GmbH
Hammweg 8–10
D-76549 Hügelsheim/Germany
Tel: +49-7229-69912-0
Fax: +49-7229-69912-20
E-Mail: info@medentika.de
www.medentika.de

Products indicated with [®] are registered brand names of the manufacturers.