


>> The new  
**flexibility** <<




**MPS**

Multi Platform Systems

TITANIUM BASE

ASC FLEX

 **MEDENTIKA**<sup>®</sup>  
A Straumann Group Brand



## ➤➤ Titanium base ASC Flex – guarantees maximum flexibility ◀◀

The titanium base for angled screw channels has been especially developed for complex prosthetics. For unfavorably positioned implants or esthetically demanding cases, it is now possible to move the screw channel in an oral direction, or to customize the chimney in four different lengths depending on the patient's needs.

Already available for the following compatible implant systems:

<b>C-Series</b>	compatible with	<b>Altatec</b>	Camlog®*
		<b>MEDENTiKA®</b>	Procone
<b>L-Series</b>	compatible with	<b>Straumann</b>	Bone Level
<b>N-Series</b>	compatible with	<b>Straumann</b>	Soft Tissue Level
<b>IPS</b>		<b>MEDENTiKA®</b>	Microcone
<b>IPS</b>		<b>MEDENTiKA®</b>	Quattrocone

The titanium base ASC Flex will be supplemented by the following compatible implant systems during the course of 2019.


<b>B-Series</b>	compatible with	<b>Bredent Medical</b>	SKY®*
<b>BS-Series</b>	compatible with	<b>BEGO Implant Systems</b>	Semados®* SC/SCX/RS/RSX/RI Semados®* S/RI
<b>CX-Series</b>	compatible with	<b>Medentis Medical</b>	ICX
<b>D-Series</b>	compatible with	<b>Altatec</b>	Conelog®*
<b>E-Series</b>	compatible with	<b>Nobel Biocare</b>	NobelReplace®* Tapered
<b>EV-Series</b>	compatible with	<b>Dentsply Implants</b>	ASTRA TECH OsseoSpeed EV
<b>F-Series</b>	compatible with	<b>Nobel Biocare</b>	NobelActive®* NobelReplace®* Conical
<b>H-Series</b>	compatible with	<b>BIOMET 3i</b>	Certain®*
<b>I-Series</b>	compatible with	<b>BIOMET 3i</b>	External Hex
<b>K-Series</b>	compatible with	<b>Nobel Biocare</b>	Brånemark System®*
<b>OT-SeriesS</b>	compatible with	<b>OSSTEM Implants</b>	TS System
		<b>HiOssen Implant®*</b>	ET System
		<b>T-Plus Implant Tech</b>	A+ Implant ST Implant
<b>R-Series</b>	compatible with	<b>Zimmer Dental</b>	Tapered Screw-Vent®*
		<b>MIS</b>	SEVEN Internal Hex
		<b>BioHorizons</b>	Tapered Internal Tapered Internal Plus Tapered Tissue Level
<b>S-Series</b>	compatible with	<b>DENTSPLY Implants</b>	ASTRA TECH OsseoSpeed®* TX
<b>T-Series</b>	compatible with	<b>DENTSPLY Implants</b>	XiVE®* S
<b>Y-Series</b>	compatible with	<b>DENTSPLY Implants</b>	ANKYLOS®* C/X

\*is a registered trademark of an independent third party

PIONEER  
TITANIUM  
BASE

One design

» ASC-Flex



MEDENTiKA®s Original «



1st Generation

**The pioneer:**

MEDENTiKA® first titanium base on the market



2nd Generation

**The evolution:**

- two different chimney heights
- adapted emergence profile




**Latest  
Generation  
ASC Flex**

**State-of-the-art:**

- angled screw channel
- 4 possible chimney heights
- optimized bonding interface
- optimized, slimmer emergence profile


compatible with all major  
implant systems

# »» MEDENTiKA® Titanium base ASC Flex ««




## IDEAL STEP WIDTH

The step width of min. 0.6 mm takes into account the requirements of a wide variety of ceramic restoration materials. This allows the safe use of press-on ceramic in accordance with the specific manufacturer's instructions.




## ROUNDED DESIGN

Reduced stress peaks, thereby protecting the ceramic restoration.




## OPTIMIZED EMERGENCE PROFILE

The optimized, even slimmer emergence profile supports and protects the soft tissue.



## INTERNAL ROTATION LOCK

The internal rotation lock receives all the material strength of the restoration, effectively avoiding predetermined breaking points. At the same time it secures the precise positioning of the hybrid abutment crown during bonding.



## BIO-PLATFORM DESIGN

Slightly tapered platform to hold the fixing material, thereby reducing the bonding gap in the gingival region.




# »» MEDENTiKA® Titanium base ASC Flex ««

## Titanium base ASC Flex


The titanium base can be used with either the angled or the straight screw channel. Thanks to the high chimney height, the titanium base offers very good support for the restoration.

### VARIABLE CHIMNEY HEIGHT




### VARIABLE CHIMNEY HEIGHT

While the chimney height of 6.5 mm also supports high restorations, it can be shortened in individual cases to 5.5/4.5/3.5 mm, reducing the vertical distance for perfect adaptation to the clinical situation.



## » Angled screw channel «




The Ball-Torx placement instrument allows angulation of up to 25°.

The maximum achievable angulation essentially depends on the respective implant bond and the abutment screw.

The angled screw channel means that the opening can be moved to areas that are less relevant in esthetic and functional respects. Accessibility can also be increased in cases where space is limited due to antagonists.

## » Snap-off point of the screw channel «



The deeper angulation point of the screw channel (in the window area) allows the point at which the screw channel emerges from the restoration to be optimized.

Series	Compatibility	Implant connection of the titanium base	Chimney height of the titanium base	Angulation of the titanium base
<b>B-Series</b>	<b>Bredent Medical / SKY®*</b>		3.5-6.5 mm	20°
<b>C-Series</b>	<b>Altatec / Camlog®*</b>	D 3.3 - D 4.3	3.5-6.5 mm	24°
		D 5.0	3.5-6.5 mm	18°
	<b>MEDENTIKA® / Procone</b>	D 3.3 - D 4.3	3.5-6.5 mm	24°
		D 5.0	3.5-6.5 mm	18°
<b>F-Series</b>	<b>Nobel Biocare / NobelActive®* / NobelReplace®* Conical</b>	RP 4./5.0	3.5-6.5 mm	20°
<b>N-Series</b>	<b>Straumann / Tissue Level</b>	NNC 3.5	3.5-6.5 mm	20°


Angulation of up to 25° is possible with all other series and implant connections.

Digital libraries are available for the following manufacturers\*:



\* to some extent this depends on the availability of the updates of the specific manufacturer.

## »» Alignment of the screw channel ««



As a rule, the screw channel is always angled in the direction of the scanbody flat (SF).

With implant systems that only allow four or fewer positioning options between abutment and implant, a second titanium base variant that is angled over the corner of the scanbody (SC) is available. This ensures greater flexibility in aligning the screw channel in the desired direction.

<b>C-Series</b>	compatible with	<b>Altatec</b>	Camlog®*
		<b>MEDENTiKA®</b>	Procone
<b>D-Series</b>	compatible with	<b>Altatec</b>	Conelog®*
<b>E-Series</b>	compatible with	<b>Nobel Biocare</b>	NobelReplace®* Tapered
<b>L-Series</b>	compatible with	<b>Straumann</b>	Bone Level
<b>N-Series</b>	compatible with	<b>Straumann</b>	Soft Tissue Level
<b>IPS</b>		<b>MEDENTiKA®</b>	Microcone
<b>IPS</b>		<b>MEDENTiKA®</b>	Quattrocone

In this case, the opening of the titanium base (angulation of the screw channel) points towards the corner of the scanbody (SC).

\*is a registered trademark of an independent third party

## SELECTION VIA THE SCANBODY:

As a rule, the screw channel is always angled in the direction of the scanbody flat (SF). With implant systems that only allow four or fewer positioning options between abutment and implant, a second titanium base variant that is angled over the corner of the scanbody (SC) is also available. This ensures greater flexibility in aligning the screw channel in the desired direction.

### Example:

#### Type SF (Scanbody Flat):

Titanium base E 1600-1-SF / Screw channel angled over the flat of the scanbody.

#### Type SC (Scanbody Corner):

Titanium base E 1600-2-SC / Screw channel angled over the right corner of the scanbody.




## SELECTION VIA THE IMPLANT CONNECTION:

An alternative to the alignment of the angled screw channel via the scanbody is alignment via the type1/type2 implant connection.

### Example:


**Type 1:** over the flat of the implant connection

**Type 2:** over the corner/cam of the implant connection



## »» Ball-Torx placement instrument ««

All titanium bases ASC Flex are screwed in with the Ball-Torx placement instrument (M 03-8 or M 10-8), guaranteeing reliable force transfer.



Contra-angled handpiece M 03-8      Manual and ratchet M 10-8



## »» MEDENTiKA® Titanium base ASC Flex ««

Product overview



# » B-Series « Bredent Medical SKY®\*

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 25 Ncm




Implant connection	NP
Chimney height	3,5-6,5
Gingiva height	0,35 mm
Article No. Type 1/SF	B 1600-1-SF
Abutment screw	B 62

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » C-Series « Altatec Camlog®\* / MEDENTIKA® Procone

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SC = Screw channel angled over the right corner of the scanbody
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 20 Ncm




Implant connection	D 3,3	D 3,8	D 4,3	D 5,0
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,35 mm	0,35 mm	0,35 mm	0,35 mm
Article No. Type 1/SC	C 1600-1-SC	C 1610-1-SC	C 1620-1-SC	C 1630-1-SC
Article No. Type 2/SF	C 1600-2-SF	C 1610-2-SF	C 1620-2-SF	C 1630-2-SF
Abutment screw	C 67	C 67	C 67	C 68

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » BS-Series « BEGO Implant Systems \*\* Semados®\* SC/SCX/RS/R SX/RI Semados®\* S/RI

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 30 Ncm



Implant connection	D 3,25/3,75-PS	D 4,1-PS	D 4,5-PS	D 5,5-PS
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,6 mm	0,45 mm	0,15 mm	0,3 mm
Article No. Type 1/SF	BS 1600-1-SF	BS 1610-1-SF	BS 1620-1-SF	BS 1630-1-SF
Abutment screw	BS 61	BS 61	BS 61	BS 61

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » CX-Series « Medentis Medical ICX

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 30 Ncm



Implant connection	D 3,75-4,8
Chimney height	3,5-6,5
Gingiva height	1,1 mm
Article No. Type 1/SF	CX 1600-1-SF
Abutment screw	CX 61

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

\*\* Products indicated with „PS“ in the implant connection are compatible with:  
BEGO Implant Systems / Semados®\* SC/SCX/RS/R SX/RI



## » D-Series « Altatec Conelog®\*


### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SC = Screw channel angled over the right corner of the scanbody
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 20 Ncm



Implant connection	D 3,3	D 3,8/4,3
Chimney height	3,5-6,5	3,5-6,5
Gingiva height	1,0 mm	1,0 mm
Article No. Type 1/SC	D 1600-1-SC	D 1610-1-SC
Article No. Type 2/SF	D 1600-2-SF	D 1610-2-SF
Abutment screw	D 62	D 62

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.



## » EV-Series « DENTSPLY Implants ASTRA TECH OsseoSpeed®\* EV


### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 25 Ncm



Implant connection	D 3,0	D 3,6	D 4,2	D 4,8	D 5,4
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	1,15 mm	1,15 mm	1,15 mm	1,15 mm	1,15 mm
Article No. Type 2/SF	EV 1600-2-SF	EV 1610-2-SF	EV 1620-2-SF	EV 1630-2-SF	EV 1640-2-SF
Abutment screw	EV 70	EV 67	EV 68	EV 69	EV 69

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.



## » E-Series « Nobel Biocare NobelReplace®\* Tapered

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Type SC = Screw channel angled over the right corner of the scanbody
- Recommended torque: 35 Ncm



Implant connection	NP 3,5	RP 4,3	WP 5,0
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,35 mm	0,35 mm	0,35 mm
Article No. Type 1/SF	E 1600-1-SF	E 1610-1-SF	E 1620-1-SF
Article No. Type 2/SC	E 1600-2-SC	E 1610-2-SC	E 1620-2-SC
Abutment screw	E 68	E 69	E 69

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.



## » F-Series « Nobel Biocare NobelActive®\* NobelReplace®\* Conical

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 15 Ncm: D 3,0  
35 Ncm: NP 3,5  
35 Ncm: RP 4,3/5,0



Implant connection	D 3,0	NP 3,5	RP 4,3/5,0
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	1,0 mm	1,0 mm	1,0 mm
Article No. Type 1/SF	F 1620-1-SF	F 1600-1-SF	F 1610-1-SF
Abutment screw	F 69	F 67	F 68

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.



## » H-Series « BIOMET 3i Certain®\*

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 20 Ncm



Implant connection	D 3,4	D 4,1	D 5,0
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,35 mm	0,35 mm	0,35 mm
Article No. Type 1/SF	H 1600-1-SF	H 1610-1-SF	H 1620-1-SF
Abutment screw	H 63	H 63	H 63

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

## » K-Series « Nobel Biocare Brånemark System®\*\*

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 35 Ncm




Implant connection	NP 3,5	RP 4,1
Chimney height	4,5-6,5	4,5-6,5
Gingiva height	0,5 mm	0,5 mm
Article No. Type 1/SF	K 1600-1-SF	K 1610-1-SF
Abutment screw	K 63	K 64

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

## » I-Series « BIOMET 3i External Hex

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 35 Ncm




Implant connection	D 3,4	D 4,1	D 5,0
Chimney height	4,5-6,5	4,5-6,5	4,5-6,5
Gingiva height	0,5 mm	0,5 mm	0,5 mm
Article No. Type 1/SF	I 1600-1-SF	I 1610-1-SF	I 1620-1-SF
Abutment screw	I 62	I 62	I 62

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

## » L-Series « Straumann Bone Level

### Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Type SC = Screw channel angled over the right corner of the scanbody
- Recommended torque: 35 Ncm



Implant connection	NC 3,3	RC 4,1/4,8
Chimney height	3,5-6,5	3,5-6,5
Gingiva height	1,0 mm	1,0 mm
Article No. Type 1/SF	L 1600-1-SF	L 1610-1-SF
Article No. Type 2/SC	L 1600-2-SC	L 1610-2-SC
Abutment screw	L 64	L 64

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » N-Series « Straumann Tissue Level

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Type SC = Screw channel angled over the right corner of the scanbody
- Recommended torque: 35 Ncm



Implant connection	NNC 3,5	RN 4,8	WN 6,5
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,6 mm	0,4 mm	0,4 mm
Article No. Type 1/SF	N 1630-1-SF	N 1610-1-SF	N 1620-1-SF
Article No. Type 2/SC	N 1630-2-SC		
Abutment screw	N 68	N 67	N 67

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » R-Series « Zimmer Dental Tapered Screw-Vent®\* MIS SEVEN Internal Hex BioHorizons Tapered Internal/Tapered Internal Plus/Tapered Tissue Level

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 30 Ncm



Implant connection	D 3,5	D 4,5	D 5,7
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,5 mm	0,4 mm	0,3 mm
Article No. Type 1/SF	R 1600-1-SF	R 1610-1-SF	R 1620-1-SF
Abutment screw	R 63	R 63	R 63

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » OT-Series « OSSTEM Implant TS-System HiOssen Implant® ET-System T-Plus Implant Tech A+ Implant, ST Implant

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 20 Ncm: M 30 Ncm: R



Implant connection	M	R
Chimney height	3,5-6,5	3,5-6,5
Gingiva height	1,0 mm	1,1 mm
Article No. Type 1/SF	OT 1600-1-SF	OT 1610-1-SF
Abutment screw	OT 62	OT 63

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# » S-Series « DENTSPLY Implants ASTRA TECH OsseoSpeed®\* TX

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Recommended torque: 15 Ncm: D 3,0 20 Ncm: D 3,5/4,0 25 Ncm: D 4,5/5,0



Implant connection	D 3,0	D 3,5/4,0	D 4,5/5,0
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	1,2 mm	1,0 mm	0,7 mm
Article No.	S 1630-SF	S 1600-SF	S 1620-SF
Abutment screw	S 70	S 71	S 69

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# T-Series DENTSPLY Implants XiVE® S

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 25 Ncm



Implant connection	D 3,4	D 3,8	D 4,5	D 5,5
Chimney height	3,5-6,5	3,5-6,5	3,5-6,5	3,5-6,5
Gingiva height	0,35 mm	0,35 mm	0,6 mm	0,6 mm
Article No. Type 1/SF	T 1600-1-SF	T 1605-1-SF	T 1610-1-SF	T 1620-1-SF
Abutment screw	T 63	T 63	T 63	T 63

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# IPS MEDENTIKA® Microcone/Quattrocone

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Type SC = Screw channel angled over the right corner of the scanbody
- Recommended torque:  
15 Ncm: NI  
25 Ncm: RI



Implant connection	NI	RI
Chimney height	3,5-6,5	3,5-6,5
Gingiva height	1,2 mm	1,2 mm
Article No. Type 1/SF	1-09-06	2-09-19
Article No. Type 2/SC	1-09-07	2-09-20
Abutment screw	1-06-04	2-06-07

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx 0-13-60 or 0-13-59. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# Y-Series DENTSPLY Implants ANKYLOS® C/X

## Titanium base ASC Flex

- angled screw channel
- Titanium Grade 5 CF
- incl. abutment screw
- Type SF = Screw channel angled over the flat of the scanbody Type
- Recommended torque: 15 Ncm




Implant connection	D 3,5-7,0
Chimney height	3,5-6,5
Gingiva height	1,2 mm
Article No. Type 1/SF	Y 1600X-1-SF
Abutment screw	

**Please note:** This is used with Scanbody 2. generation. To screw in the titanium base ASC Flex you need the Placement instrument Ball-Torx M 03-8 or M 10-8. To select the desired direction of the angled screw channel, please consider the Instruction for use.

# Tools

## Placement instrument Ball Torx

- Hardened stainless steel



Version	Contra-angle	Manual and ratchet
Type		
Article No.	M 03-8	M 10-8



Publisher: MEDENTiKA® GmbH  
Hammweg 8-10  
76549 Hügelsheim  
Tel: +49 (0)7229 69912-10  
info@medentika.de  
www.medentika.com

Design: The WeberFink GbR  
graphic design studio  
www.weberfink.de

Date: February 2019

We are certified:  
**DIN EN ISO 13485**  
**Medizinprodukterichtlinie 93/42/EWG,**  
**Annex II**

**CE0483**

Technical changes and errors reserved.

You can find the Instructions for use and warranty conditions on the website [www.medentika.com](http://www.medentika.com)  
More information on the warranty can also be requested directly from the manufacturer



MEDENTiKA® GmbH  
Hammweg 8-10  
76549 Hügelsheim  
info@medentika.de  
www.medentika.com

»» Passion  
for **Precision** ««