



WILLI GELLER
Creation
SIMPLY BRILLIANT

Creapearl Denture Teeth





Type A Anatomic - Upper anterior with natural rootshaped emergence profile

The denture tooth is a perfect anatomical copy of a natural tooth.

Type V Vario - Upper anterior teeth for individual contouring of denture teeth

The tooth length and appearance can be influenced with or without grinding of the neck by designing the gum. The single tooth can get an individual appearance.

Type U - Lower anterior

The same can be applied to the lower anterior teeth which can be arranged in a similar way like the type V Vario. The user is able to determine the length of the tooth, and there is no need for a large tooth stock.



WILLI GELLER *Creapearl*

THE BRILLIANCE OF CREAPEARL:

- Anterior denture tooth to satisfy the highest demands
- Dynamic depth of colour thanks to highly chromatic colouring of the dentine
- Highly opalescent incisors guarantee a natural look and can be customised without colour loss
- The naturally fluorescent and translucent material gives a brilliant light transmission
- Special designed interapproximal areas for a natural red-white look
- Expressive, natural-looking shapes with deep interdental fastening panel to counter shadowing on individual arrangements!
- Manual production process with an individual note



WILLI GELLER
Creation
 SIMPLY BRILLIANT

C R E A P E A R L Combination chart		Colours
Anterior upper	Anterior lower	
A71	U3	A1
A73	U7	A2
A75	U7	A3
A78	U7	A3,5
A79	U5	B2
V61	U1	B3
V63	U3	C3
V65	U5	
V68	U9	
V69	U5	

Distributor
 Creation Willi Geller International GmbH
 Koblacherstraße 3, 6812 Meiningen, Austria
 Tel. +43 (0)5522 76784, Fax +43 (0)5522 76784-59
 info@creation-willigeller.com, www.creation-willigeller.com

Manufacturer
 KLEMA Zahnproduktion GmbH
 Koblacherstr. 3a, 6812 Meiningen, Austria



YOUR LOCAL DEALER